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The screenshot shows the 'アラーム設定' (Alarm Settings) tab. It is divided into three color-coded sections: yellow for bed status, cyan for breathing, and red for pulse. Each section has input fields for a threshold value and a timing (minutes and seconds) for when an alarm is triggered.

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The screenshot shows the 'メール内容設定' (Email Content Settings) tab. It features six text input fields arranged in a 2x3 grid. The top row is for '離床' (bed leaving), '呼吸(しきい値以上)' (breathing above threshold), and '脈(しきい値以上)' (pulse above threshold). The bottom row is for '着床' (bed getting), '呼吸(しきい値以下)' (breathing below threshold), and '脈(しきい値以下)' (pulse below threshold). Each field is currently empty.

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
The screenshot shows the 'メール設定' (Email Settings) tab. It is divided into two main sections: 'サーバー設定' (Server Settings) and 'メールアドレス設定' (Email Address Settings). The server settings section includes a port number field, an 'SSL使用' (Use SSL) checkbox, and fields for SMTP server, user name, and password. The email address settings section includes fields for the sending email address, and three recipient email addresses (TO, CC, CC). A 'メール送信テスト' (Test Email) button is located at the bottom right of the email address section.

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
The screenshot shows the 'アラーム設定' (Alarm Settings) tab. It is divided into three color-coded sections: yellow for '離床着床' (Bed Exit/Bed Entry), cyan for '呼吸' (Respiration), and red for '脈拍' (Pulse). Each section contains input fields for 'しきい値' (Threshold) and '分' (Minutes) and '秒' (Seconds) for timing. All values are currently set to 0. Buttons for '登録' (Register) and '終了' (End) are at the bottom.

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The screenshot shows the 'メール内容設定' (Email Content Settings) tab. It features six text input fields arranged in a 2x3 grid. The top row is for '離床' (Bed Exit), '呼吸(しきい値以上)' (Respiration above threshold), and '脈(しきい値以上)' (Pulse above threshold). The bottom row is for '着床' (Bed Entry), '呼吸(しきい値以下)' (Respiration below threshold), and '脈(しきい値以下)' (Pulse below threshold). Buttons for '登録' (Register) and '終了' (End) are at the bottom.

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The screenshot shows the 'メール設定' (Email Settings) tab. It is divided into two main sections: 'サーバー設定' (Server Settings) and 'メールアドレス設定' (Email Address Settings). The server settings section includes a 'ポート番号' (Port Number) field, an 'SSL使用' (Use SSL) checkbox, and fields for 'SMTPサーバー' (SMTP Server), 'ユーザー名' (Username), and 'パスワード' (Password). The email address settings section includes fields for '送信メールアドレス' (Sending Email Address), '受信メールアドレス(TO:)' (Receiving Email Address (TO:)), and '受信メールアドレス(CC:)' (Receiving Email Address (CC:)). A 'メール送信テスト' (Email Send Test) button is located at the bottom right of the email address section. Buttons for '登録' (Register) and '終了' (End) are at the bottom.

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